



Session 1:
EEGLAB
(in MATLAB)

Daniel Bennett
UniSA, March 2016






Introduction
& set-up of computers



EEGLAB
Program overview



EEGLAB
A whistle-stop tour



Introduction
& set-up of computers



EEGLAB
Program overview



EEGLAB
A whistle-stop tour




Introductions (briefly)

Daniel Bennett

PhD Candidate, University of Melbourne

I research decision making under uncertainty, using a combination of neuroimaging methods (EEG, fMRI) and computational cognitive modelling.

Obligatory embarrassing EEG headshot:





Introduction to EEGLAB

EEGLAB is an open-source toolbox for preprocessing and analysis of EEG data.

MATLAB-based

Freely available (though MATLAB is not)

This 1-hour session will cover the basics of EEGLAB.


Features and limitations

Pointers toward some useful references

No MATLAB knowledge required.

We will largely use the GUI, not the command line

Information about command-line access to data will also be provided for those so inclined.



Setting Up Computers

Download the eeglab folder from this link: bit.ly/1WVrClt

Unzip it into your My Documents > MATLAB directory

In MATLAB, click 'Set Path' > 'Add folder and subfolders'

If prompted to save a different path file ('pathdef.m'), save this to your MATLAB folder too.




Introduction
& set-up of computers



EEGLAB
Program overview



EEGLAB
A whistle-stop tour



Introduction to EEGLAB

EEGLAB is an open-source toolbox for preprocessing and analysis of EEG data.

MATLAB-based

Freely available (though MATLAB is not)

This 1-hour session will cover the basics of EEGLAB.


Features and limitations

Pointers toward some useful references

No MATLAB knowledge required.

We will largely use the GUI, not the command line

Information about command-line access to data will also be provided for those so inclined.



History of EEGLAB

1997-2001: Initially developed as an Independent Components Analysis toolbox at Salk Institute, La Jolla, USA.


2001-present: Maintained by the Schwartz Centre for Computational Neuroscience, UCSD.

Website: <http://sccn.ucsd.edu/eeglab/>

Citation: Delorme, A., & Makeig, S. (2004). EEGLAB: an open source toolbox for analysis of single-trial EEG dynamics including independent component analysis. *Journal of Neuroscience Methods*, 134(1), 9-21.

Version:

- Current version is v.13 (MATLAB 2014a and later)
- Old version is v. 12 (pre-2014a)



Pros and Cons

Pro:

- Access to single-trial data
- Open-source and non-proprietary code
- Flexible batch scripting
- Extensive documentation
- Surprisingly many inbuilt features
- Ability to create new modules
- Data visualisation

Cons:

- Memory requirements
- Unlovely GUI
- Code stability
- User support
- Data visualisation

For us, the most important considerations were to have control over every step of the analysis, and to have theory-neutral and flexible access to data.



Introduction
& set-up of computers



EEGLAB
Program overview




EEGLAB
A whistle-stop tour



Getting started

1. Open MATLAB
2. Type `eeglab` at the command line
3. Press enter.






Getting started

These warnings are distracting, and kind of ugly...

```
Command Window
>> eeglab
eeglab: options file is C:\Users\dbennett1\eeg_options.m
EEGLAB warning: there can be only one EEGLAB window, closing old one
EEGLAB: adding "dipfit" v2.3 (see >> help eegplugin_dipfit)
EEGLAB: adding "firfilt" v1.6.1 (see >> help eegplugin_firfilt)
Warning: The file 'C:\Users\dbennett1\TOOLBOXES\eeglab13_5_4b\eeglab.m' could not be cleared because it contains MATLAB code that is
currently executing.
> In eeglab>updatemenu (line 1416)
   In eeglab (line 368)
   In eeglab (line 1064)
Warning: The file 'C:\Users\dbennett1\TOOLBOXES\eeglab13_5_4b\eeglab.m' could not be cleared because it contains MATLAB code that is
currently executing.
> In eeglab>updatemenu (line 1416)
   In eeglab (line 368)
   In eeglab (line 1064)
Warning: The file 'C:\Users\dbennett1\TOOLBOXES\eeglab13_5_4b\eeglab.m' could not be cleared because it contains MATLAB code that is
currently executing.
> In eeglab>updatemenu (line 1416)
   In eeglab (line 368)
   In eeglab (line 1064)
You are using the latest version of EEGLAB.
fx >>
```

1. Enter the following at the command line: `edit eeglab`
2. Insert a % symbol at the beginning of line 1398: `%clear functions;`
3. Resave the file.




Importing Data

1. File > Load existing dataset > 'eeglab_data.set'.

NB: The file format for existing eeglab datasets is **.set**. If we had wanted to import directly from the raw datafile, we could have used File > Import data > Using EEGLAB functions and plugins.

If you look in the data folder, you will notice that there are two datafiles called eeglab_data. The .set file contains metadata, and the .fdt file contains the eeg recordings.




Visualising Data

Plot > Channel data (scroll)

Things to note:

1. Timestamp
2. Vertical scale
3. Events
4. Stacking/unstacking
5. Rejecting data




Edit > Channel locations...

Visualising Data

Plot > Channel spectra and maps

Plotting frequency range: [2 64]




Visualising Data

Tools > Filter the data > Basic FIR filter

Lower edge = 58

Higher edge = 62

‘Notch filter the data...’ = ticked




Visualising Data

Tools > Filter the data > Basic FIR filter

Lower edge = 58

Higher edge = 62

‘Notch filter the data...’ = ticked




Visualising Data

Tools > Filter the data > Basic FIR filter

Lower edge = 58

Higher edge = 62

‘Notch filter the data...’ = unticked





Accessing Data

Data can be easily accessed from the command line, in the EEG structure:

```
>> EEG  
>> EEG.data
```


Continuous data like in this file are stored in an $A \times B$ matrix, where
A is the number of electrodes,
B is the number of data points

```
>> EEG.data(5, 1:128)  
  
>> plot(EEG.data(5, 1:128))  
  
>> EEG.data(5, :) = 0
```


Analysing Data


We will now open up a file which has already been 'epoched' (divided into segments for analysis)

1. File > Clear dataset(s)
2. File > Load existing dataset > 'eeglab_data_epochs_ica.set'.
3. Plot > Channel data (scroll)



Note that data is now discrete, not continuous (try clicking on a segment).

Note also that time 0 in a segment always corresponds to a green line.



Analysing Data

Unlike continuous data, epoched data are stored in an $A \times B \times C$ matrix, where

A is the number of electrodes,

B is the number of data points per epoch

C is the number of epochs

As a result, we can access single-trial data from the command line as follows:


```
>> EEG.data(:, :, 1)
```

And, if we so desire, we can access information about all the events that happen within an epoch as follows:

```
>> EEG.epoch(1)
```


Plotting Averaged Data


1. Plot > Channel ERPs > With scalp maps



Plotting Averaged Data

1. Plot > Channel ERPs > In scalp/rect array

Click on individual channels to enlarge




Plotting Averaged Data

1. Plot > Channel ERP image

Channel = 27

Time limits (ms) = -500 1500




Plotting Averaged Data

1. Plot > Channel ERP image

Channel = 27

Time limits (ms) = -500 1500

Smoothing = 0



Plotting Averaged Data

1. Plot > Channel ERP image


Channel = 27

Time limits (ms) = -500 1500

Epoch-sorting field = 'latency'

Event type(s) = 'rt'

Smoothing = 10





Some Final Comments

Moving from point-and-click to batch scripting:

```
>> eegh
```

EEGLAB is tremendously useful for some things, especially for the purposes of MVPA.

It has some limitations in other respects, though.
For instance, no function to detect amplitude of a peak.

In many cases, these limitations are addressed by separate plugins

- ERPLAB (Lopez-Calderon & Luck, 2014)
- sLORETA

Useful resource: http://sccn.ucsd.edu/wiki/EEGLAB_TUTORIAL_OUTLINE